



Equipment

- One Collection Bottle
- One Paediatric urine collection bag

Instructions

1. Wash the genital area thoroughly. Dry thoroughly.
2. Remove the adhesive protector from the plastic bag.
3. Carefully stick the adhesive tab over the genital area covering the infant's penis [boys] or labia (skin folds) [girls], carefully avoiding the anus.
4. After the urine has been passed, carefully remove the bag and pour the urine into the urine container.
5. Sample should be delivered as soon as possible on the day of collection. If there is a delay in delivery, please refrigerate.